

Geometric accuracy of 3D-printed vs. milled dental crowns: A surface-partitioned analysis of occlusal, intaglio and marginal regions

Giulia Pascoletti¹ Francesco Bianconi¹ Giuliana Baiamonte² Michele Cali²
Nicola Scotti³ Elisabetta M. Zanetti¹

¹Department of Engineering, University of Perugia, 06125, Perugia, Italy

²Electric, Electronics and Computer Engineering Department, University of Catania, Catania, 95125, Italy

³Department of Surgical Sciences, Dental School of the University of Torino, 10126 Turin, Italy

Overview

1. Background and objectives
2. Materials
3. Methods
4. Results
5. Discussion
6. Conclusions and future work

1. Background and objectives

2. Materials

3. Methods

4. Results

5. Discussion


6. Conclusions and future work

Dental crowns: what they are why they are used

Custom-made caps that permanently cover a damaged, weak or decayed tooth

Functions:

- **Restoration** (to restore a damaged or worn-down tooth)
- **Protection** (to protect a weak or partially broken tooth)
- **Cosmetics** (to cover a misshapen tooth or an implant)



Source: Wikimedia Commons
(CC BY-SA 3.0)

Dental crowns: how they are made (1/2)

Traditional (analogue) methods

- Lost-wax casting
- Heat pressing



Subtractive methods

- CAD/CAM milling



Additive manufacturing (3D printing)

- Digital light processing (DLP)
- Lithography-based ceramic manuf. technology (LCM)
- Stereolithography (SLT/TSLA)



Dental crowns: how they are made (2/2)

Subtractive manufacturing

- High geometric accuracy
- Good mechanical strength



Additive manufacturing

- Less material waste compared to milling
- Intricate geometries possible



Dental crowns: functional regions

Occlusal region (blue):


- Chewing (external) surface of the crown
- Touches the opposing teeth when biting

Intaglio region (purple):

- Inside (internal) surface of the crown
- Adheres to the prepared tooth

Marginal region (white):

- Boundary between the intaglio and occlusal
- Guarantees sealing



Dental crowns: why precision matters

For the **occlusal** region:

- Maintains bite and jaw alignment
- Improves masticatory efficiency

For the **intaglio** region:

- Determines fit and retention

For the **marginal** region:

- Prevents micro-leakage and bacteria infiltration

Objectives

- To investigate **geometric accuracy** (trueness) of **milled** vs. **additive** manufactured dental crowns
 - Global accuracy analysis
 - Separate analysis for the occlusal, intaglio and marginal regions

1. Background and objectives

2. Materials

3. Methods

4. Results

5. Discussion

6. Conclusions and future work

Reference geometry–virtual model

Physical starting sample:

- Intact upper molar, extracted for periodontal reasons

Acquisition:



- Intraoral 3D scanner Trios 3 (3Shape)

Generation of the crown model:

- Dental system CAD software (3Shape)
- Uniform thickness offset (1.5 mm)

Virtual sample generated:

- Triangulated STL model (34374 faces)



Physical samples

Six groups of specimens:

- Milled samples (two groups, M1–M2)
- 3D-printed samples (four groups, A1–A4)

Group	Material	Technique	Num. samples
M1	Composite	Milling	7
M2	Ceramic	Milling	8
A1	Composite	AM/DLP	7
A2	Ceramic/composite	AM/TSLA	6
A3	Ceramic	AM/LCM	9
A4	Ceramic	AM/LCM	16

Milled samples

Group M1

- Material: BRILLIANT Crios (Coltene Group)
- Milling machine: inLab MC XL (Dentsply Sirona)



Group M2

- Material: KATANA Zirconia STML (Kuraray Noritake Dental Inc.)
- Milling machine: inLab MC XL (Dentsply Sirona)



3D-printed samples (1/2)

Group A1

- Material: VarseoSmile Crown plus (BEGO GmbH & Co.)
- Printer: SprintRay Pro S (SprintRay Inc.)
- Thickness: 50 μm



Group A2

- Material: Irix Max (DFAB)
- Printer: DFAB Desktop 3D (DWS Systems)
- Thickness: 50 μm



3D-printed samples (2/2)

Group A3

- Material: 5Y-TZP zirconia (Tosoh Corporation)
- Printer: LCM (Admatec Admaflex 130)
- Thickness: 50 μm



Group A4

- Material: 8Y-TZP zirconia (Tosoh Corporation)
- Printer: LCM (Admatec Admaflex 130)
- Thickness: 30 μm



1. Background and objectives

2. Materials

3. Methods

4. Results

5. Discussion

6. Conclusions and future work

Accuracy estimation procedure


1. Sample digitalisation
2. Shape registration
3. Surface segmentation
4. Evaluation of the goodness-of-fit
5. Statistical analysis

Sample digitalisation

3Shape Trios 3 system (same used for generating the reference geometry)

Declared specs:

- Weight: 340 g
- Capture rate: 1875 images/sec
- Precision: $4.5 \pm 0.9 \mu\text{m}$
- Trueness: $6.9 \pm 0.9 \mu\text{m}$





Source: Hokayem et al., 2025
(CC BY 4.0)

Shape registration

Performed with **CloudCompare** (Open Source Project)

- a. Preliminary manual alignment (translation and rotation)
- b. Fine registration by iterative closest point (ICP)
 - Stopping criterion: $< 15 \mu\text{m}$ diff. between two consecutive iterations



Evaluation of the goodness-of-fit (trueness)

Local performance metric (node level):

- **Signed distance** δ between each node of the sample's mesh and the closest triangle in the reference geometry
 - **Positive** for deviations in the same direction as the surface normal (angle between 0° and 90°)
 - **Negative** for deviations in direction opposite to the surface normal (angle between 90° and 180°)

Global performance metric (surface level):


- Root mean square (**RMS**) of δ

Segments analysed:

- **Whole surface**
- **Occlusal**, **intaglio** and **marginal** regions separately


Surface segmentation: identification of the marginal region

- a. Manual generation of an approximating closed polyline (CloudCompare)
- b. Polyline over-sampling (Matlab-custom routine)
- c. Identification of the closest nodes on the mesh (Matlab-custom routine)



Surface segmentation: occlusal and intaglio regions

- Occlusal: all nodes not belonging to the marginal line and with the normal vectors pointing **outward**
- Intaglio: all nodes not belonging to the marginal line and with the normal vectors pointing **inward**



1. Background and objectives

2. Materials

3. Methods

4. Results

5. Discussion

6. Conclusions and future work


Trueness

Group	Surface segment			
	Overall	Intaglio	Occlusal	Marginal
M1	79.4 $\pm$ 4.74	47.0 $\pm$ 8.97	95.7 $\pm$ 9.33	41.2 $\pm$ 5.76
M2	94.3 $\pm$ 63.3	125.6 $\pm$ 105.9	55.5 $\pm$ 19.2	83.5 $\pm$ 45.3
A1	112.0 $\pm$ 12.5	156.4 $\pm$ 23.8	72.9 $\pm$ 16.1	83.2 $\pm$ 19.4
A2	78.2 $\pm$ 25.6	100.6 $\pm$ 44.1	56.3 $\pm$ 14.1	74.5 $\pm$ 35.7
A3	111.6 $\pm$ 11.1	127.8 $\pm$ 13.7	98.3 $\pm$ 18.7	110.6 $\pm$ 27.0
A4	118.1 $\pm$ 13.2	121.6 $\pm$ 22.9	115.1 $\pm$ 13.0	114.9 $\pm$ 25.1

- RMS values for δ pooled by group (mean $\pm$ std., μm)


Deviation colourmaps: milled samples

- Intaglio surface:
 - General material **deficiency**
- Occlusal surface:
 - Fairly contained overall deviation
 - Concentrated material **excess** (M1)



Deviation colourmaps: 3D-printed samples

- Intaglio surface:
 - General material **excess**
- Occlusal surface:
 - Material **deficiency** in large areas of A3–A4
 - No clear trend for A1–A2



1. Background and objectives

2. Materials

3. Methods

4. Results

5. Discussion

6. Conclusions and future work

Clinical fitness of the samples

- Cut-off values for trueness:
 - 100–150 μm generally considered acceptable
 - 50–100 μm generally considered good
- Results aligns well with the current literature
- Nearly all groups (except A1) achieved acceptable values for all the surface segments
- Group M1 attained good values for all the surface segments

Comparing methods

- All 3D-printed groups resulted in **material deficiency** at the intaglio surface
- A3-A4 yielded **material excess** on the occlusal surface
- M1 generated localised material notches at the occlusal surface

1. Background and objectives

2. Materials

3. Methods

4. Results

5. Discussion

6. Conclusions and future work

Conclusions

- The accuracy of 3D-printed crowns was clinically acceptable on the whole
- Milling on composite (M1) was good for the intaglio and marginal regions, however milling on ceramics (M2) performed better at the occlusal region
- The intaglio surface was more critical for AM than the occlusal region

Future work

- Further validate the results on **larger datasets**
 - More reference teeth
 - Larger number of physical samples
- Investigate reasons behind **material excess/deficiency** in specific zones

Aknowledgements

Partially supported by **Fondazione Cassa di Risparmio di Perugia**:

- Call Research 2025 (Innovative research projects, group A)
- Project **MODUS-PRO**–*Automated and personalised design of prostheses and ortheses* (Ref. 2025.0582)



Thank you for your attention

Any questions?